Recurrent Fever **SYPHILIS** 

Frambæsia Sleeping Sickness

# GALYL

# of GALYL were given in one of our Service Hospitals.

"No case has given the slightest cause for anxiety, and the clinical results are very satisfactory."—Staff-Surgeon's Report.

#### IMPROVED INTRAMUSCULAR INJECTION.

The GALYL is in the form of a Sodium Salt in a weak solution of Glucose. This is much easier to use and more rapidly absorbed than the oily emulsion. It is in every way identical with the intravenous, — efficacy, dose, price, etc., — only differing in administration.

Thorough and exhaustive tests in various parts of the Kingdom have proved it to be equal in every way with the intravenous; moreover, patients may resume and continue their occupations without inconvenience after injection of intramuscular GALYL.

·Simplified Wassermann Test.

# LATAPIE

#### Method of the Pasteur Institute.

The Latapie reaction for Syphilis which only requires three test tubes, has marked advantages over other methods. It simplifies the technique and makes errors almost impossible.

This statement is based upon extensive practical experience embracing more than 100,000 examinations conducted at the Pasteur Institute (since 1908) with blood collected at the Institute itself or sent in from outside sources.

#### Accurate—Reliable—Uncomplicated.

The Antigen is supplied in tubes, which is diluted 22-26 times (according to indications on each tube).

It keeps good and fully active for several years if kept at room temperature and protected from light.

Full and Complete Instructions with each Tube.

LITERATURE CLINICAL REPORTS, and PRICE LIST on REQUEST.

THE ANGLO-FRENCH DRUG CO., Ltd., Gamage Building, Holborn, London, E.C. 1.

(late M. Bresillon & Co.).

Telegrams: "AMPSALVAS, LONDON."

#### HYPERTHERMAL THERAPY.

Thone: Holborn 1311.

### ambrine

Ambrine is the only Authorised preparation of Dr. Barthe de Sandfort. We are compelled to state this, as there are various imitations.



May be applied with camel-hair brush, but our specially designed spray greatly facilitates the application of Ambrine, especially in the case of wounds or extensive surfaces.

Ambrine forms a dressing which INSTANTLY alleviates pain—promotes rapid healing—forms a sure protection against iffection of the wound—can be removed without pain, hamorrhage, or injury to the newly-formed tissues. The healing leaves no scars or contractions.

The most efficient means and method of ensuring active Hyperthermal-Therapy to the practical treatment of Rheumatism, Gout, Arthritis, Yaricose Ulcers, and similar affections.

## URASEPTINE

## The Most Powerful and Effective Urinary Antiseptic.

URASEPTINE is a granulated product entirely soluble in water, its bases being Piperazine, Urotropine, Helmitol, Benzoates of Sodium and Lithium. It contains 60 centigrams (10 grs.) of active matter to each teaspoonful.

Dose. - 2 to 6 teaspoonfuls daily.

It purifies the Urine, and this action is due to its three principal properties:

- I. It is a URINARY ANTISEPTIC.
- 2. A SOLVENT of URIC ACID and of PHOSPHATES.
- 3. A MILD NON-TOXIC DIURETIC.

INDICATIONS. — Arthritism, Gout, Gravel, Hepatic and Renal Colic, Rheumatism, Calculus, &c., Phosphaturia, Urinary Antisepsis, Pyelitis, Bacteriuria, Cystitis, Prostatitis, Urethritis, Pyuria, Urinary Abscess, Vesical Catarrh, &c.