

TWO CASES THAT MAY MAKE YOUR NEXT ONE EASIER



Emergency, on-the-scene radiography and screening is precisely the situation in which the Victor Model F-3-35 finds its greatest usefulness. Two well-balanced carrying cases contain a complete, practical, portable-mobile X-ray unit consisting of the well-known Model F-3 X-ray Unit and the Model 35 Tube Stand.

Enjoy the many benefits of real X-ray portability—benefits as close to you as

the nearest plug point of suitable voltage and frequency. And this, at no loss to your consulting room or hospital auxiliary service, since the F-3-35 is also a fine-quality mobile unit. Yet, when the need arises, remove the F-3 from its tube stand, easily and quickly dismantle and pack the few component parts snugly in their twin carrying cases, and your X-ray is ready for "Flying Squad" service.

Still, the F-3-35 is more than merely a portable-mobile X-ray unit. It is easy to use; safe and silent in operation. Its powerful, wide-range effectiveness is amply demonstrated by the short exposures recommended in the unit's comprehensive technique charts for virtually all regions of the body. Moreover, the busy F-3-35 is moderately priced for so versatile a unit and is available on a convenient plan adjusted to today's budgets.

A free, factual booklet makes clear the many practical, exclusive features of the F-3-35. For a copy, fill up the convenient coupon and post it, today.

---WITHOUT OBLIGATION---

VICTOR X-RAY CORPORATION, Ltd.

15-19 CAVENDISH PLACE, LONDON, W. 1

Please send a copy of the free booklet on the Victor F-3-35 Portable X-ray Unit.

Name.....

Address.....

Town.....

KG-2


Wheat-starch conversion by **BENGER'S FOOD**




Benger's Food is a specially prepared wheaten powder containing adequate quantities of natural pancreatic enzymes of which the proteolytic enzyme known as trypsin and the starch converting amylases are of first importance.

The wheaten powder serves two highly important functions (1) it acts as conveyor of these enzymes (2) in preparation of the Food by the admixture of hot milk it is converted into dextrins, dextri-maltoses and maltose.

When the hot milk meets the cold mixture the resulting temperature is correct for the enzyme action. Within two minutes the starch of the wheaten powder has been made soluble so that prepared Benger's Food contains



only 0.5% of insoluble starch. As the digestion proceeds the trypsin modifies the milk protein so that when it comes in contact with the gastric juices, it separates into fine flocculi, as compared to the heavy tough curd formed by milk itself. Both actions take place whilst the Food is too hot to drink.



For Infants, Invalids and those suffering from severe digestive weakness, the time for pre-digesting Benger's

Food should be extended to 15 minutes and upwards to 45 minutes. At the end of this period the wheaten powder is converted to dextrins, dextri-maltoses and maltose, and the milk proteins modified so that the finest possible curd is formed on contact with the gastric juices.

Sole Manufacturers: **BENGER'S FOOD LTD.**

Holmes Chapel, Cheshire, Eng.

The Benger Laboratories are always at the disposal of the medical profession in any of their dietetic problems. Physicians' samples of Benger's Food are always available.