

WICKHAM HEAD LAMP

This popular and effective Head Lamp is now standard equipment used by the Armed Forces. Synthetic non-perishable cable, chrome plated finish, and sound construction give an instrument of lasting use.

Plugs into standard Gowlland Battery Handles.

Catalogue No. 19.

Gowlland
electric diagnostic instruments

Made
in
England

Obtainable from all Surgical Supply Houses.

for RAPID DISINFECTION ARGYROL Silver Vitellin plus NATURAL EPHEDRINE IMMEDIATE DECONGESTION

in the treatment of

**COLDS • ACUTE OR CHRONIC RHINITIS • CATARRH
SINUSITIS • ASTHMA • HAY FEVER • CONJUNCTIVITIS**

The combination into a stable solution of Silver Vitellin and natural laevorotatory Ephedrine—a powerful and lasting congestion-relieving agent—in normal saline makes Argotone the ideal treatment for nasal affections and their complications. ARGOTONE is an ethical preparation not advertised to the Public, and can be prescribed with every confidence.




Medical samples available on request

ARGOTONE

NASAL & EYE DROPS

A stabilised compound of Argylol 1% Ephedrine 0.9%
in Normal Saline

RONA LABORATORIES LIMITED, 159, FINCHLEY ROAD, LONDON, N.W.3



SCIENTIFIC GLASSWARE

**provides long
and efficient
service.**

PYREX Brand Scientific Glassware is supplied only through Laboratory Furnishers, but illustrated catalogue and two free copies of our Chemist's Notebook will be sent direct on application to us.

Ask for PYREX Brand and see that you get it!
JAMES A. JOBLING & CO. LTD.
Wear Glass Works, SUNDERLAND.

T.47.E.

THE protection afforded by the use of PYREX Brand Scientific Glassware in all routine dispensing, and in all meticulous research work, is truly remarkable.

Not only is this Glassware immune from the effects of sudden thermal changes, but it is also highly resistant to all acids (except hydrofluoric and glacial phosphoric).

By virtue of the amazingly low co-efficient of expansion of .000032, the structure of PYREX Brand Scientific Glassware can be made more robust than that of ordinary glass.

This gives it additional mechanical strength, which enables it to resist the physical shocks of everyday usage . . . thereby saving much of the expense of glassware replacements.

MOIST HEAT IN THE TREATMENT OF BOILS AND CARBUNCLES

The moist heat of an ANTIPHLOGISTINE pack plus its hygroscopic and osmotic action is of definite value in the treatment of boils and carbuncles.

Applied comfortably hot the moist heat of Antiphlogistine is maintained for many hours. Pain and soreness are lessened and the local reparative processes are stimulated.

ANTIPHLOGISTINE may be used with chemo-therapy.

FORMULA:

Glycerinum.	45.000%	Methylum Ess ...	00.002%
Iodum	00.010%	Eucalypti Ess ...	00.002%
Acid Boricum ...	00.100%	Kaolin	54.864%
Acid Salicylicum	00.020%		100.000%
Menthae Ess	00.002%		

Antiphlogistine
TRADE MARK

The Denver Chemical Mfg. Co.
LONDON, N.W.9