



Aseptic

Surgery

LONDON



HOSPITAL

# 'ULTRATAN' CATGUT

(Morley's Process)


'Ultratan' Catgut gives absolute sterility in ligatures. It is prepared under the strictest aseptic precautions in the special laboratories which have been set up for the purpose at the London Hospital and is in use in the 15 operating theatres of the establishment. The unique feature is that sterilisation is effected immediately after the splitting of the intestine before the gut is twisted. This ensures that no contaminating debris, which invariably escapes sterilisation, is enclosed within the turns, later to become the focus of infection at the site of operation.

**In tubes for  
Immediate  
Surgical Use**

—  
In  
10, 20, 30 & 40 day.  
Standard thicknesses.

Nos. 00, 0, 1, 2,  
3, 4 & 5

—  
Three tubes in a box.



**Sterility** of every batch is verified before issue both by Aerobic and Anaerobic Tests.

**Absorption** is practically a standard factor and can be guaranteed in every grade.

**Tensile strength** is such that even in the finest gut there is a wide margin of safety.

Descriptive literature on request to  
the sole distributors:

**Allen & Hanburys Ltd.**  
48 Wigmore Street, W. 1.

**In 100 feet  
lengths for  
Hospitals**

—  
10, 20, 30 & 40  
day gut.

In boxes of 100 ft.

Nos. 00, 0, 1, 2,  
3, 4 & 5.

—  
Requires External  
Sterilisation.

