

# BULLOCK AND REYNOLDS' E C L E C T I C   I N H A L E R ,

AS RECOMMENDED BY DR. MORELL MACKENZIE,

*Physician to the Hospital for Diseases of the Throat; and Assistant Physician to the London Hospital.*

This Apparatus can be used for all Medicated Vapours, and may be employed in all positions of the patient. It requires *no effort in inspiration*, insures the *thorough medication* of the Vapour, and accurately *maintains the temperature* desired.

	WHOLESALE.	RETAIL.
Price with Thermometer - - - -	7s. 6d.	10s. 0d. each
Ditto with Stand and Spirit Lamp complete - -	9s. 6d.	12s. 6d. ,,

DESCRIPTIVE HANDBILLS ON APPLICATION.

**S. MAW, SON, & THOMPSON,**

10, 11, & 12, ALDERSGATE STREET, LONDON, E.C.

## MAW'S DOUBLE-VALVED EARTHENWARE INHALER,

For the inhalation of the Vapour of Hot Water, either alone or impregnated with Ether, Chloroform, Henbane, Creasote, Vinegar, etc., in affections of the Throat and Bronchial Tubes, Asthma, Consumption, etc.

This apparatus is especially adapted for Chloroform, Ether, and other active and volatile substances, as from its construction the exact quantity required may be inhaled. The *lower* valve opens only during *inhalation*, whilst the *exhaled* breath passes freely from the *upper* valve, instead of through the vessel, dissipating a large portion of the medicated vapour, as is the case in the use of most other inhalers.

Price to the Trade, 4s. 6d. each; Retail, 7s.

DESCRIPTIVE HANDBILLS ON APPLICATION.

**S. MAW, SON, & THOMPSON,**

10, 11, & 12, ALDERSGATE STREET, LONDON, E.C.

## DR. NELSON'S IMPROVED EARTHENWARE INHALER

Will be found a most efficient Apparatus for the inhalation of the Vapour of Hot Water, either alone or impregnated with Ether, Chloroform, Henbane, Creasote, Vinegar, etc., in affections of the Throat and Bronchial Tubes, Asthma, Consumption, etc.

Price to the Trade, 3s 6d each; Retail, 5s. 6d.

DESCRIPTIVE HANDBILLS ON APPLICATION.

**S. MAW, SON, & THOMPSON,**

10, 11, & 12, ALDERSGATE STREET, LONDON, E.C.